

PRODUCT INFORMATION

PRODUCT NAME: Large-Insulated-Medical-Kit--
 6xTS-Plastic-Container
PRODUCT CODE: STP-4G57-EPS67-IP2-6xTS
 (IS-IP-X6)

PRODUCT TYPE: Infectious Substances and
 Biological Substances

DESCRIPTION: '620/650 Insulated Specimen Transport Pack' with fibreboard outer box, EPS insulating box, EPS cage support system with lid and 6 x 'Hazpak' tube sets (Primary Receptacle and Secondary Package) for transport with chilling mediums or at ambient temperature. IATA COMPLIANT.

This package can be used with 1 or up to 6 x tube sets as per the client's requirements.

Quantity of chilling medium will depend on the number of tube sets used.

NOTE: If transporting UN 2814 or UN 2900 place the 6.2 Class label over the UN 3373 label, and mark proper shipping name over 'Biological Substance Category B' marking.

PHOTOGRAPHS



Packaging System

FALCON TEST ENGINEERS REPORT No: 0594

CONFIGURATION OF PACKAGING:

'Hazpak' tube sets are inserted into the EPS cage support system and the EPS insulating box closed. The top flaps of the fibreboard outer box are closed and sealed using self adhesive tape. Add the required amount of absorbent material if specified. Label the box as required. If using dry ice, the box should be pierced in the position marked to allow for venting.

MARKINGS & APPROVAL:

Approval No.: 3378
UN Markings: UN 4G/CLASS 6.2/###AUS/3378/HAZPAK
Packaging Group: N/A
S.G.: N/A
Internal Test Pressure: 180 kPa at Ambient and -95 kPa at +55°C & -40°C

PRODUCT DIMENSIONS:

Dimensions: 380 (L) x 380 (W) x 385 (H) mm **Mass of Empty Package:** 2197 grams

Nominal Volume:	0.06 m ³ per unit
Bulk Equivalent:	14.3 kg cubic (@ 250 kg/m ³)

LIST OF COMPONENTS:

Fibreboard outer box	Class 9.1 label
EPS insulating box	Documents enclosed envelope
Hazpak tube sets	Class 6.2 label
EPS cage support system with lid	UN 3373 label

ASSEMBLY OF PACKAGING:

Assemble corrugated fibreboard outer box and secure bottom flaps using self-adhesive tape. Place the EPS insulating box into the fibreboard outer box, and remove the lid of the EPS insulating box. Using the jig affix the EPS cage support system base to the centre of the inside of the EPS insulating box base, and the EPS cage support system lid to the centre of the inside of the EPS insulating box lid. Close the EPS insulating box. Mark the top of the EPS insulating box lid and the inside of the flap of the fibreboard box next to the mark. **Do Not** place the tube sets in to the EPS cage as we supply bagged tube sets separately. Place the required labels on top of the EPS insulating box. Close and seal the top flaps of the fibreboard box with self-adhesive tape.

OUTER PACKAGE:

Product Code:	4G57	Qty:	1
Description:	Corrugated fibreboard box.		
Internal Dimensions:	375 (L) x 375 (W) x 395 (H) mm		

INSULATING BOX:

Product Code:	EPS67	Qty:	1
Description:	EPS insulating box.		
Dimensions:	370 (L) x 370 (W) x 375 (H) mm		

PRIMARY RECEPTICLE AND SECONDARY PACKAGE:

Product Code:	RF-IP2-TS	Qty:	6
Description:	Hazpak tube set.		
Dimensions:	42 (D) x 175 (H) mm	Internal Diameter:	33 mm

INSERT:

Product Code:	EPS-CAGE-LID	Qty:	1
Description:	EPS cage support system with lid.		
Dimensions:	180 (L) x 130 (W) x 235 (H)		

LABELS:

Product Code:	LAB-6.2-100x100	Qty:	1
Description:	100 x 100 mm Class 6.2 label.		<input type="text"/>
Product Code:	LAB-UN3373-100x100	Qty:	1
Description:	100 x 100 mm UN 3373 label.		<input type="text"/>
Product Code:	LAB-9.1-100x100	Qty:	1
Description:	100 x 100 mm Class 9 Miscellaneous label.		<input type="text"/>
Product Code:	Documents Envelope	Qty:	1
Description:	Documents enclosed envelope.		<input type="text"/>



Components