

## PRODUCT INFORMATION

**PRODUCT NAME:** MINI MEDICAL KIT      **PRODUCT TYPE:** Infectious Substances and  
**PRODUCT CODE:** STP-4G4-IP2-1xTS      Biological Substances  
**DESCRIPTION:** '620/650 Ambient Temperature Specimen Transport Pack' with fibreboard outer  
 box, fibreboard stabilising insert, 2 x EPS stabilising inserts and 'Hazpak' tube set  
 (primary receptacle and secondary package) for transport at ambient temperature.  
 IATA COMPLIANT

**NOTE:** If transporting UN 2814 or UN 2900 place the 6.2 Class label over the UN 3373 label, and mark proper shipping name over 'Biological Substance Category B' marking.

### PHOTOGRAPH



### Components

**FALCON TEST ENGINEERS REPORT No.:** 8512A

### CONFIGURATION OF PACKAGING:

'Hazpak' tube set is inserted into the fibreboard stabilising insert. The top flaps of the fibreboard stabilising insert are folded over and the EPS stabilising insert placed on top. Fibreboard outer box is closed and secured by using the self-locking tab. Add the required amount of absorbent material if specified. Label the box as required.

**(Make sure that the labels or any other markings are not obliterated or obscured in any way).**

Note: The use of other packaging methods or components may render the approval of this package invalid

### MARKINGS & APPROVAL:

**Approval No.:** 3506  
**UN Markings:** UN 4G/CLASS 6.2/\*\*/AUS/HAZPAK3506  
**Packaging Group:** N/A  
**S.G.:** N/A  
**Internal Test Pressure:** 180 kPa at Ambient and -95 kPa at +55° C and -40° C

### PRODUCT DIMENSIONS:

**Dimensions:** 145 (L) x 130 (W) x 220 (H) mm      **Bulk Equivalent:** 1.1 kg (@ 250 kg/m<sup>3</sup>)  
**Nominal Volume:** 0.004 m<sup>3</sup> per unit      **Mass of Empty Package:** 254 g