

HAZPAK PTY LTD 44 Stephen Rd, Dandenong South Vic. 3175 Australia Postal Address: P.O. Box 4000, Dandenong South Vic. 3164 Australia Tel: (03) 9706 8058 Fax: (03) 9706 7593 Intl. Tel: 61 3 9706 8058 E-mail: hazpak@hazpak.com.au A.B.N. 67 054 007 162

## **PRODUCT INFORMATION**

PRODUCT NAME: Small-Medical-Kit-2L-Plastic- PRODUCT TYPE: Biological Substances Pouch-

PRODUCT CODE: STP-4G6-IP5-1x2L

**DESCRIPTION: '650 Ambient Temperature Specimen Transport Pack'** with fibreboard outer box, re-sealable plastic bag (Primary Receptacle), **'Hazpak'** Medi-Pouch (Secondary Package) and 4 x absorbent pads for transport at ambient temperature. IATA COMPLIANT.

**Note:** The absorbent material supplied absorb up to 200 mL of product. Additional absorbent pads can be purchased separately if required.

## PHOTOGRAPHS



Packaging System

FALCON TEST ENGINEERS REPORT No.:

4023

## CONFIGURATION OF PACKAGING:

Four absorbent pads and the re-sealable plastic bag are placed inside the '**Hazpak**' Medi-pouch. The pouch is sealed following printed instructions and rolled lengthways and placed inside fibreboard box. The top flap of the fibreboard outer box is closed and secured by using the self-locking tab. Add the required amount of absorbent material if specified. Label the box as required.

(Make sure that the labels or any other markings are not obliterated or obscured in any way).

Note: The use of other packaging methods or components may render the approval of this package invalid

MARKINGS & APPROVAL:			
Approval No.:	N/A		
UN Markings:	N/A		
Packaging Group:	N/A		
S.G.:	N/A		
Internal Test Pressure:	-95 kPa		
PRODUCT DIMENSIONS:			
Dimensions:	120 (L) x 120 (W) x 240 (H) mm	Bulk Equivalent:	0.9 kg (@ 250 kg/m <sup>3</sup> )
Nominal Volume:	0.004 m <sup>3</sup> per unit	Mass of Empty Package:	115 g