

HAZPAK PTY LTD 44 Stephen Rd, Dandenong South Vic. 3175 Australia Postal Address: P.O. Box 4000, Dandenong South Vic. 3164 Australia Tel: (03) 9706 8058 Fax: (03) 9706 7593 Intl. Tel: 61 3 9706 8058 E-mail: hazpak@hazpak.com.au A.B.N. 67 054 007 162

PRODUCT INFORMATION

PRODUCT NAME:	Small-Medical-Kit-30mL-Test-	PRODUCT TYPE:	Infectious Substances and
	Tube		Biological Substances
PRODUCT CODE:	STP-4G7-IP2-6xTS		

DESCRIPTION: '620/650 Ambient Temperature Specimen Transport Pack' with fibreboard outer box, EPS cage support system with lid and 6 x 'Hazpak' tube sets (Primary Receptacle and Secondary Package) for transport at ambient temperature. IATA COMPLIANT NOTE: If transporting UN 2814 or UN 2900 place the 6.2 Class label over the UN

3373 label, and mark proper shipping name over 'Biological Substance Category B' marking.

PHOTOGRAPHS



Packaging System

FALCON TEST ENGINEERS REPORT No.: 8507

CONFIGURATION OF PACKAGING:

'Hazpak' tube sets are inserted into the EPS cage support system and the EPS cage support system lid is placed on top. Close and seal the box with 48 mm pressure sensitive self adhesive tape in the same manner that the base of the box has been sealed. Add the required amount of absorbent material if specified. Label the box as required.

(Make sure that the labels or any other markings are not obliterated or obscured in any way)

Note: The use of other packaging methods or components may render the approval of this package invalid.

MARKINGS & APPROVAL:				
Approval No.:	3377			
UN Markings:	UN 4G/CLASS 6.2/##AUS/HAZPAK 3377			
Packaging Group:	N/A			
S.G.:	N/A			
Internal Test Pressure:	180 kPa at Ambient and -95 kPa at +55°C & -40°C			
PRODUCT DIMENSIONS:				
Dimensions:	190 (L) x 135 (W) x 270 (H) mm	Bulk Equivalent:	1.7 kg (@ 250 kg/m ³)	
Nominal Volume:	0.007m ³ per unit	Mass of Empty Package:	928 grams	