



A.B.N. 67 054 007 162  
HAZPAK PTY LTD 44 Stephen Rd, Dandenong South Vic. 3175 Australia  
Postal Address: P.O. Box 4000, Dandenong South Vic. 3164 Australia  
Tel: (03) 9706 8058 Fax: (03) 9706 7593 Intl. Tel: 61 3 9706 8058  
E-mail: hazpak@hazpak.com.au

## PRODUCT INFORMATION

**PRODUCT CODE:** RF-IP2-2L-95-IP5-1L (DS-IS-BIOJAR2)      **PRODUCT TYPE:** Infectious Substances and Biological Substances Refurbishment

**DESCRIPTION:** **Refurbishment Only - Primary and Secondary Packaging System** with re-sealable plastic bag (Primary Receptacle), 2 L HDPE round jar with 95 mm screw closure and o-ring gasket (Secondary Package) and 4 x absorbent pads approved for use with 'Hazpak' 650 Ambient Temperature and Insulated Specimen Transport Packaging systems only.

### PHOTOGRAPHS



**Secondary Package**



**Packaging System**

**FALCON TEST ENGINEERS REPORT No:** 3340 and 3489

### CONFIGURATION OF PACKAGING:

Four absorbent crystal pads are placed inside the 'Hazpak' 2 L HDPE jar and the re-sealable plastic bag is placed on top of the absorbent crystal pads. The jar is closed and sealed by screwing close the lid.

### MARKINGS & APPROVAL:

**Approval No.:** 3033  
**UN Markings:** N/A  
**Packaging Group:** N/A  
**S.G.:** N/A  
**Internal Test Pressure:** 180 kPa at Ambient and -95 kPa at +55°C and -40°C

### PRODUCT DIMENSIONS:

**Dimensions:** 230 (H) x 119 (Dia) mm      **Mass of Empty Package:** 275g  
**Nominal Volume:** 0.003 m<sup>3</sup> per unit  
**Bulk Equivalent:** 0.8 kg (@ 250 kg/m<sup>3</sup>)  
**Internal Orifice Dimensions:** 85 mm (Dia)