

A.B.N. 67 054 007 162

HAZPAK PTY LTD 44 Stephen Rd, Dandenong South Vic. 3175 Australia Postal Address: P.O. Box 4000, Dandenong South Vic. 3164 Australia Tel: (03) 9706 8058 Fax: (03) 9706 7593 Intl. Tel: 61 3 9706 8058

E-mail: hazpak@hazpak.com.au

PRODUCT INFORMATION

PRODUCT CODE: STP-4G57-EPS67-IP2-3x2L- **PRODUCT TYPE**: **Biological Substances**

95

DESCRIPTION: '650 Insulated Specimen Transport Pack' with fibreboard outer box, EPS insulating

box, 5 x EPS stabilising inserts, 3 x re-sealable plastic bags (Primary Receptacle), 3 x

'Hazpak' 2 L HDPE round jars with 95 mm screw closure and o-ring gasket

(Secondary Package) and 12 x absorbent pads for transport with chilling mediums or

at ambient temperature. IATA COMPLIANT

This package can be used with 1 or up to 3 x 2L jars as per the client's requirements.

Quantity of chilling medium will depend on the number of 2L jars used.

NOTE: The absorbent material supplied absorb up to 200 mL of product. Additional

absorbent pads can be purchased separately if required.

PHOTOGRAPHS



Secondary Packaging



Packaging System

FALCON TEST ENGINEERS REPORT No: 4040

CONFIGURATION OF PACKAGING:

Four absorbent pads and the re-sealable plastic bag are placed inside each of the 'Hazpak' 2 L HDPE jars. The 2 L HDPE jars are inserted into the EPS insulating box and closed. The top flaps of the fibreboard outer box are closed and sealed using self adhesive tape. The box is labelled as required. If using dry ice, the box should be pierced in the position marked to allow for venting.

MARKINGS & APPROVAL:

Approval No.: 3378

UN Markings: 4G/CLASS 6.2/##/AUS/3378/HAZPAK

Packaging Group: N/A S.G.: N/A

180 kPa at Ambient and -95 kPa at +55° C and -40° C **Internal Test Pressure:**

PRODUCT DIMENSIONS:

380 (L) x 380 (W) x 395 (H) mm Mass of Empty Package: **Dimensions:** 2535 grams

> 0.057 m³ per unit **Nominal Volume**